



Together through life sciences.

pCLAMP workshop

Jeffrey Tang, PhD 2013

#### **Axon Conventional Electrophysiology Family**


- pCLAMP software
  - Clampex-data acquisition
  - Clampfit-Data analysis
  - Axoscope-Data continuous monitoring
- Digitizer
  - Digidata 1440A-Analog-to-Digital converter
- Patch-clamp amplifiers
  - Axopatch 200B
  - MultiClamp 700B
  - AxoClamp 900A



#### Together through life sciences.

©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.

00000











#### Flow of recording signal



#### **Electrophysiology rig**



#### Together through life sciences.



#### **Measurement Techniques**

- Extracellular voltage recording
  - Single unit and field potential recording in brain or brain slices, electrocardiagrams, encephalograms and oculograms
- Voltage clamp
  - Voltage is held constant and the current passing through the cell membrane is measured
  - Patch-clamp recording
- Two electrode voltage clamp
  - Voltage is held constant through one electrodes and current is measured with the other electrode
- Current clamp
  - Current is held constant and the corresponding membrane voltage of the cell is measured
- Intracellular sharp electrode recording
  - Measurement of action potentials
  - Bridge balance for the pipette resistance
- Discontinuous clamp
  - Amplifier divides its time between passing current and recording voltage
- Ion-selective electrodes and electrochemistry
  - Measurement of small changes in ion, neurotransmitter her hand hormone concentrations in tissues and cells

©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.

#### **MolDev Axon CEP Products**



#### **The position of Molecular Devices**

- Market leader
  - Axon brand
    - over 25 years
  - High quality
    - ultra low noise amplifiers
    - High speed, low noise digitizers
  - Reliable
    - Lasting in a good condition for many years
  - Technical support
    - Professional assistance

- Large install base (17,000)
  - Academic research labs
  - Pharmaceutical/biotechnology drug discovery labs
  - >10K amplifiers sold
  - >10K digitizers sold
  - >10K copies of pCLAMP



#### **Customer support activities**

- Scientific conferences
  - Society for Neuroscience
  - Biophysical Society
  - FENS
- Online webinars
  - Getting the Most Out of pCLAMP series

- pCLAMP Workshops at university
  - Duke University
- University loaner program/Support training courses
  - Cold Spring Harbor marine biology lab
  - Woods Hole marine biology lab
  - Australian Course in Advanced Neuroscience (ACAN)
  - Patch clamp workshop in Singapore
  - More...



#### Together through life sciences.





# A Walkthrough of Protocol Editor in Clampex



Molecular Devices

# **Golden Tips that help**

- Study manual
- Study manual again
- Ask someone who uses pCLAMP
- Clampex tutorial
- Online HELP
- Knowledge base
- Webinar tutorials
- Technical Support



#### Together through life sciences.

# Agenda

Edit Protocol - (untitled)	
Mode/Rate Inputs Outputs Trigger Statistics Comments Math Waveform Stimulu	is

- Acquisition Mode
- Inputs
- Outputs
- Trigger
- Statistics
- Comments
- Math
- Waveform





# **Feature Highlights**

- Acquisition modes
- Scope trigger
- Holding level overrides
- Digital OUT holding pattern
- Stimulus file
- Digital bit pattern
- Pre-sweep train
- P/N leak subtraction
- User list
- Membrane test between sweeps




#### Together through life sciences.

# Telegraphs

 Clampex can receive and incorporate arrange of "telegraphed" amplifier settings as the variable gain, lowpass filter, whole-cell capacitance compensation.



# Telegraph setting in AxoPatch<sup>™</sup> 200B



#### Together through life sciences.



# **Telegraph setting in AxoPatch<sup>™</sup> 200B**

Lab Bench 🔀	Lab Bench 🛛 🗙
Input Signals         Digitizer Channels       Signals         Ipatch       Add         Inalog IN #0       Image IN #0         Analog IN #1       Image IN #0         Analog IN #3       Image IN #0         Analog IN #3       Image IN #0         Analog IN #4       Image IN #0         Analog IN #5       Image IN #0         Analog IN #5       Image IN #0         Analog IN #6       Image IN #0         IN 0       Rename         IN 0       Scaling         Signal units:       p A         Default       Scale factor (V/pA):         0.001       Scale Factor         Offset (pA @ 0 V):       0         Software RC Filter       Hardware Signal Conditioning         Lowpass (kHz):       Image Information         Highpass (Hz):       Image Information         Frequency (Hz):       10000	Input Signals       Digitizer Channels       Signals         Image IN #0       Add         Analog IN #1       Add         Analog IN #1       Image IN         Analog IN #2       Image IN         Analog IN #3       Image IN         Analog IN #4       Image IN         Analog IN #4       Image IN         Analog IN #5       Image IN         IN 0       Rename         IN 0_KK       Image IN         Signal units:       Image IN         Scaling       Signal units:         Signal units:       Image IN         IN 0       Scale Factor         Offset (pA @ 0 V):       Image IN         Software RC Filter       Hardware Signal Conditioning         Lowpass (Hz):       Image IN         Gain:       5         Cmap IN       Image IN         Gain:       5         Cmap IN       Image IN         Image IN       Image IN         Image IN       <
OK Cancel Help	OK Cancel Help
Together through life	SCALED OUTPUT OUTPUT GAIN ( ( ) OVLD x 10 x 20 x 2 x 100 x 2 0 x 1 x 200 x 1 x 200 x 2 0 x 1 x 200 x 200 x 1 x 200 x 200

Molecular Devices

## **Telegraph setting in MultiClamp<sup>™</sup> 700B**

Telegraphed Instrument	🔩 MultiClamp 700B: ( 00103420 ) 📃 🗖 🔀	Lab Bench 🛛 🗙					
Input Channels       Telegraphed Instrument         Analog IN #0       Analog IN #1         Analog IN #1       Analog IN #2         Analog IN #2       Analog IN #2         Analog IN #2       Analog IN #3         Analog IN #3       Analog IN #4         Analog IN #5       Figure 1         Analog IN #5       Figure 1         Analog IN #6       Figure 1         Analog IN #7       Figure 1         Analog IN #6       Figure 1         Analog IN #6       Figure 1         Analog IN #7       Figure 1         Analog IN #8       Figure 1         Analog IN #8       Figure 1         Analog IN #9       Figure 1         Analog IN #9       Figure 1         Analog IN #10       Continue trial when mode, scale factor or units change.	Channel 1       I (pA)         -2.6       -6.4         Resistance       Ims         Mode       VC         VC       IC         Ext       V-Clamp 1         V-Clamp 1       V-Clamp 2	Input Signals       Output Signals         Digitizer Channels       Signals         Analog IN #0       Add         Analog IN #1       Imprime         Analog IN #2       Imprime         Analog IN #3       Delete         Analog IN #4       IN 0         Analog IN #5       IN 0_KK         In 0_KK       Ipatch					
Telegraph Connections         Choose the digitizer channels to which telegraphs are connected:         Gain:       Telegraph Input 0 I Cm:         Frequency:       Telegraph Input 1 I	Holding:       0 mV         Seal Test:       100 Hz         Cp Fast:       3.850 pF         Image: Comparison of the state of t	IN 0       Scaling       Signal units:       p ▼ A       Default       Scale factor (V/pA):       0.001       Scale Factor (V/pA):       0					
Amplifier output: Primary output 1  Serial number: 00103420  Scan	Auto       Correction: 0 % Prediction: 0 %         Disable if oscillation detected         Primary Output: Membrane Current (2.5 V/nA)         Gain: 5       Bessel: 10 kHz       AC: DC       Scope: Bypass         Output Zero       Lock Subtraction	Software RC Filter       Hardware Signal Conditioning         Lowpass (kHz):       5         Highpass (Hz):       1         Telegraphs					
Output Channels Choose the Analog OUT channel to which the command signal is connected: Command: Analog OUT #0 OKCancelHelp	0 mV         Auto         10.00 MΩ         Auto           Secondary Output: Membrane Potential (10 mV/mV)         Gain: 1         Lowpass Filter: 10 kHz           Pulse         10 mV         10 ms         Zap         500 μs         Rf: 500 MΩ	Gain:         5         Cm (pF):         0.000           Frequency (Hz):         10000         Scale factor (V/nA):         0.5           OK         Cancel         Help					




#### Together through life sciences.

# Lab Bench setting in AxoPatch<sup>™</sup> 200B


🛜 Clampex - [Scope	- Demo	Digitizer]	
📑 File Edit Acquire	· View	Configure Tools Window Help	
📸 🗣 🛃 🔈		Digitizer CyberAmp •	
Episodic Crnd 0 (mV)	3 C <sub>4</sub>	Membrane Test Setup	
·80     ◆       Cmd 1 (mV)       ·80       A0 #2 (mV)       0       ▲0 #3 (mV)       0       Digital OUTs       3       7       4       Sampling (Hz)       10000	5	Sequencing Keys  Program Options Lab Book Options Toolbars File Registration  Voice Tags	Lab Bench       X       Lab Bench         Input Signals       Output Signals       Input Signals       Input Signals         Digitizer Channels       Signals       Digitizer Channels       Signals         Input Signals (N #1)       Imprime       Analog IN #2       Analog IN #3         Analog IN #2       Analog IN #3       Delete       Imprime         Innemb       Rename       Ipatch       Rename         IN 0       Rename       Ipatch       Rename         IN 0       Rename       Cmd 0       OUT 0         Scaling       Signal units:       p < A
(100 μs) IN 1 Lowpass (kHz)	-5- - -10 10		Scale factor (V/pA): Offset (pA @ 0 V): 0 Voltage Clamp ne signals for each digitizer channel Inits and scale factors
		2. 10 set t	inits and scale factors

Molecular Devices

#### Together through life sciences.



# Lab Bench setting in AxoPatch<sup>™</sup> 200B



#### Together through life sciences.



# Lab Bench setting is telegraphed in MultiClamp<sup>™</sup> 700B

Telegraphed Instrument	Lab Bench 🔀
Input Channels       Telegraphed Instrument         Analog IN #1       Analog IN #2         Analog IN #2       Analog IN #2         Analog IN #3       Amalog IN #4         Analog IN #4       Analog IN #4         Analog IN #5       Analog IN #5         Analog IN #6       Analog IN #6         Analog IN #7       Analog IN #8         Analog IN #8       Analog IN #8         Analog IN #8       For simultaneously switch protocols, configure special Sequencing Keys (V-Clamp or I-Clamp IN).         Continue trial when mode, scale factor or units change.	Input Signals Digitizer Channels Analog IN #0 Analog IN #1 Analog IN #2 Analog IN #3 Analog IN #5 Analog IN #5 Analog IN #6
Telegraph Connections         Choose the digitizer channels to which telegraphs are connected:         Gain:       Telegraph Input 0 •         Frequency:       Telegraph Input 1 •	IN 0 Signal units: p ▼ A Default Scale factor (V/pA): 0.001 Scale Factor Assistant
Amplifier Confiduration Amplifier output: Primary output 1 Serial number: 00103420 Scan	Software RC Filter Lowpass (kHz): 5 Highpass (Hz): 1 CyberAmp
Output Channels Choose the Analog OUT channel to which the command signal is connected: Command: Analog OUT #0 -	Telegraphs Glain: 1 Cm (pF): 0.000 Frequency (Hz): 10000 Scale factor (V/nA): 0.5
OK Cancel Help	OK Cancel Help




#### Together through life sciences.

# **Protocol Editor**

🛜 Clampex - [	Scope - Demo Digitize	r]													
📑 File Edit	Acquire View Configu	Edit Protocol - (unt	itled)								×				
	New Protocol	Mode/Rate Inputs	Outputs Tri	Edit Protocol - (untitled	I)										X
	Open Protocol	-Acquisition Mode-		Mode/Rate   Inputs   Outp		ioner Ì S	tatistics	Comme	nte Í Mat	h Wa	veform	Stimulus	1		
n car car a la	Save Protocol	C Gap-free		modernate   inputs   out	uts j n	igger   o	tausues	Comme				ounaias	1		-1
	Save Protocol As	C Variable-le	enath events	Waveform Analog OUT:	Cmd 0		Info								
23 G_k	Edit Protocol Waveform Preview			Analog Waveform -					🔲 Digit	al Outpul	ts —				
	Export protocol ( pCL4	Trial Hierarchy	0	🖲 Epochs 🔿 Stimu	ilus file				- <u>R</u> 4	Active hig	gh logic fa	or digital t	rains	Info	
Episodic Cmd 0 (mV)		Trial delay (s):		Intersweep holding le		Jse holdir	ng 🔻	7    F	Inter	awaan b	it pattern:		e holding		
-80	Record	Runs/trial:	1	Intersweep holding le	ver. ju	536 110101	ig		Inter	sweep n	it pattern.	103	s noiding		
Cmd 1 (mV)	Re-Record	Sweeps/run:	10	Epoch Description	A	B	C	D	E	F	G	Н		J	
-80	View Only			Туре	Step	Off	Off	Off	Off	Off	Off	Off	Off	Off	
A0 #2 (mV)	Pause-View	Sweep duration (s):		Sample rate	Fast	Fast	Fast	Fast	Fast	Fast	Fast	Fast	Fast	Fast	
0	Pause	First holding: 3.1 ms	Epochs: 193.8 ms	First level (mV) Delta level (mV)	112 -20	0	0	0	0	0	0	0	0	0	
AO #3 (mV)	Stop Repeat	31 samples	1938 san	First duration (ms)	100	0	0	0	0	0	0	0	0	0	
0	Write Last	Sampling Rate per S	Signal	Delta duration (ms)	0	0	0	0	0	0	0	0	0	0	
Digital OUTs	Autotrigger			Digital bit pattern (#3-0)	1111	0000	0000	0000	0000	0000	0000	0000	0000	0000	
ЗГГГГ		Fast rate (Hz):	10000	Digital bit pattern (#7-4)	0000	0000	0000	0000	0000	0000	0000	0000	0000	0000	
7 🗆 🗆 🗆	Time Tag	Slow rate (Hz):	10000	Train rate (Hz) Pulse width (ms)	0	0	0	0	0	U	0	U	U	0	
Sampling (Hz)	Comment Tag	j 510W Tate (H2).	1.0000	· · · ·	0	U	U	0 0	-		0	0	0	0	
10000 🚔	Voice Tag	Space available is 6	00153 sweep	Number of sweeps = 10				Alloca	ated time:		11	)6.2 of 20	JU ms		
(100 µs)	1 Conditioning.pro	Allow automatic a		Stimulus File											
– IN 0 Lowpass (kHz)_	2 Baseline.pro	Niow automatic t	andiyala in oli k	Summary											
				Juninary											
Highpass (Hz)			1	Channel #0 Channel #1	Chan	nel #2 _ (	Channel ‡	13		ternate V	Vaveform	s 🔲 Alt	ernate D	igital Outpu	uts
	-10000	OK (	Cancel												
<u> </u>	10		ito u	OK Cance	el	He	elp		quisition n Episodic		on		Upda	te Preview	1
		otocol Ed	ILOF:						Episodic	sumulau	011				-
-	lis a	the "ce	ntraľ	" for proto	col		itir								_
								·9·			_				
The second base of the	ugh life sciences.											_		lecula	

#### Together through life sciences.

## **Acquisition Mode?**





Molecular Devices

# Mode/Rate

Edit Protocol - (untitled)							
Mode/Rate	Inputs Outputs Trigger Statistics Comments Math Waveform Stimulus						

- Passive Acquisition
  - Gap-free
  - Fixed-length
  - Variable-length
  - High-speed oscilloscope

- Stimulus Acquisition
  - Episodic stimulation



# **Gap-free Mode**

- Data are acquired continuously
- No gaps in data file
- single-channel or minis recordings




 Together
 through
 life
 sciences.

 ©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.

# **Fixed-length Events Mode**


- Constant data segment above the threshold and pre-trigger portion are recorded
- Action-potential spikes or other constantwidth events recordings







#### **Fixed-length Event Mode**








# **Variable-Length Events Mode**


- Data above the threshold, and for pre- and post-trigger periods are recorded.
- Variable recorded segment
- "Bursting" data recording





#### Together through life sciences.

#### **Variable-Length Events Mode**




#### Together through life sciences.



# **High-Speed Oscilloscope Mode**

- Like an oscilloscope
- Data above the threshold, and for pre- and post-trigger periods are recorded.





#### ©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.

Together through life sciences

# **Passive acquisition modes**

dit Protocol - (untitled)
Mode/Rate Inputs Outputs Acquisition Mode © Gap-free © Variable-length eve Check this option to acquire data until the disk is full, or until acquisition is stopped manually
Trial Length         O Use available disk space         O Duration (hh:mm:ss):       00:00:05         = 0.10 MB         Space available = 57:03:16 (hh:mm:ss) (28493 MB).
Sampling Sampling rate per signal (Hz): Sweep length per event per signal (Hz): Sweep
Specify the sampling rate of the analog-to-digital conversion
Allow automatic analysis in other programs Enable Configure / Automatic Analysis in Clampfit to DK Cancel Help

#### Together through life sciences.




# **High-speed Oscilloscope Mode**

Edit Protocol - (untitled) 🔀	
Mode/Rate Inputs Outputs Trigger Statistics Comments Math Waveform Stimulus	
C Gap-free C Fixed-length events C Episodic stimulation C Variable-length events C High-speed oscilloscope	
Trial Length Averaging Options	
C Use available disk spa_Type of Average	
Duration (sweeps):     Cumulative	
Space available = 102698 C Most recent Weighting (%): 10 🚔 ~ 10 sweeps in average	
Sampling 🔽 Undo File	
Sampling rate per signal   Update every (sweeps): 1, 2, 5	
Sweep length per event C Always prompt to allow undo	
Prompt on premature termination only	
Averaging OK Cancel Help Default	
Allow automatic analysis in other program Check this option to save the raw average sweep OK Cancel He	sweeps and









- Provide a command waveform and record responses, in fixed-length sweeps
- Analog waveform, holding level and/or digital pulses are outputs
- Special features include pre-sweep trains, online leak current subtraction, online event detection and statistics, and an on-line derived-math channel.









Molecular Devices

#### Together through life sciences.




#### Together through life sciences.

#### Episodic Stimulation Mode---Holding Periods



#### Together through life sciences.





#### Together through life sciences.



#### Inputs

🖉 Clampex - [				•			
Edit Ad			Window Help	 	1 1	 ( (	1
🚔 🗣 🗳	De 🖪		09	<u> </u>	ŧ 🌌 -		
	<b>6</b>	844	€ 12-12				
<u>୧</u> ୦ ଜ <sub>ା</sub> ଜ <sub>ା</sub>	2 C <sub>3</sub> C <sub>4</sub>	ତ <sub>୍ର</sub> ାତ <sub>୍ର</sub> ାତ <sub>୍ମ</sub> ାତ୍	8 G <sub>A1</sub> G <sub>A2</sub> G	Gal Cal Ca2	EA3		
Episodic Cmd 0 (mV)							
	<u> </u>						
Cmd 1 (mV) 0 €	€ € 20000 € -20000						
Cmd 2 (nA)	-20000 3456780 100000 100000 10000 100000 10000 10000 10000 10	-					
Cmd 3 (mV) 0 €							
Digital OUTs	중 은 -10000						
	₹ 2 0 -10	-					
Sampling (Hz) 10000 🚔	S ( 10 ≥ 0 ≤ 10	-					
(100 μs)	-10						
IN 0 Lowpass (kHz)	≅ 2 0 <u>≤ -10</u>						
□ 5 🔮   Highpass (Hz)		-					
	= <u>-10</u> ∞ ( 10 ∞ ( 10 ∞ ( 10) ∞ (						
- Telegraphs —— Filter: 10 kHz	-10						
Gain:1 Cm:0pF	≊ 2 0 ≝ 2 -10	-					
Ra: 0M							
Elapsed time		-					
0:00:00 Run Sweep		-					
		-					
	v 0 10						
	<u>∑</u> ≥ 0 -10						
		o I	(ms)	50			1

#### Together through life sciences.



# Q & A

- Q: Why the selected input name in the Lab Bench does not show up in the scope window?
- A: You need to select the input name in the Input tab of the Edit Protocol.




### **Outputs**

Ð	Edit Protocol	- (untitled)				Þ		
N	Mode/Rate Inputs Outputs Trigger Statistics Comments Math Waveform Stimulus							
Analog OUT Channels								
	Channel #0:	Cmd 0	-	Range (mV): -2	200.00 to 199.99 at 20.0 mV/V	Voltage-clamp		
Channel #1: Cmd 1			-	Range (mV): -2	Range (mV): -200.00 to 199.99 at 20.0 mV/V			
Channel #2: Cmd 2			-	Range (nA): -10.00 to 10.00 at 1.0 nA/V Current		Current-clamp		
	Channel #3:	Cmd 3	•	Range (mV): -2	200.00 to 199.99 at 20.0 mV/V			
I I	-Analog OUT Ho	olding Levels			Check to specify	one digital output		
	Cmd 0 (mV):			-70		· · ·		
Cmd 1 (mV):		0	channel to go high (and stay high) during the entire length of a trial					
	Cmd 2 (nA):			0				
	Cmd 3 (mV):			0				
	- Digital OUT Hol	lding Pattern		765	4 3 2 1 0			
l								






## **Digital Out Holding Pattern**





#### Together through life sciences.



#### Together through life sciences.



## **Holding Level Overrides**

🗟 Lab Bench 🛛 🗙	🗈 Lab Bench 🛛 🔀 🕯
Input Signals Output Signals	Input Signals Output Signals
Digitizer Channels Signals	Digitizer Channels Signals
Analog OUT #0 Cmd 0 Add	Analog OUT #0 Add
Analog OUT #1 Cmd 0 Analog OUT #2 L_clamp Delete	Analog OUT #2 I_clamp Delete
Analog OUT #3 OUT 0 V_clamp Rename	V_clamp Rename
VC Cmd	
Cmd 0	Cmd 0
Signal units: m 💌 V Default	Signal units: m 💌 V Default
Scale factor (mV/V): 20 Scale Factor	Scale factor (mV/V): 20 Scale Factor
Holding level (mV):	
Analog holding level is disabled because of the setting in the Configure / Overrides dialog.	To set the analog holding level in the protocol, use Configure / Override.
Digital OUT Channels	Digital OUT Channels
7 6 5 4 3 2 1 0 Holding pattern:	7 6 5 4 3 2 1 0
Digital holding pattern is disabled because of the setting in the Configure / Overrides dialog.	To set the digital holding level in the protocol, use Configure / Override.
Set digital OUT bit high during acquisition: 🔲 Digital Bit: 4 💌	Set digital OUT bit high during acquisition: Digital Bit: 4
OK Cancel Help	OK Cancel Help



#### Together through life sciences.

## **Holding Level Overrides**

Clampex - [Scope - Demo	o Digitizer]	
📓 File Edit Acquire View (	Configure <mark>Tools Window He</mark>	elt (i
<b>F B B</b>	Digitizer CyberAmp 🛛 🕨	
	Lab Bench Telegraphed Instrument	
<u>୧</u> ୁ ାଟ <sub>ା</sub> ଟ <sub>ୁ</sub> ାଟ <sub>ୁ</sub>   ଟ <sub>ୁ</sub>   ଟ <sub>ୁ</sub>   ଟ୍ରୁ   ଟ୍ରୁ	Overrides	🗈 Overrides 🛛 🔀
		<ul> <li>Check to override individual protocol settings:</li> <li>Use analog holding levels from Lab Bench</li> <li>Use digital holding pattern from Lab Bench</li> <li>Use most recent file comment</li> <li>Use this sampling rate (Hz): 1000</li> <li>Keep existing Scope window size and position</li> <li>OK</li> <li>Cancel</li> <li>Help</li> </ul>



#### Together through life sciences.

# Q & A

 Q: How do I change the command voltage for more than 200 mV

Analog OUT Cha	og OUT Channels	
Channel #0:	Cmd 0	■ Range (mV): -200.00 to 199.99 at 20.0 mV/V

 A: You need to change the scale factor and external sensitivity of analog output signal on the Lab Bench



### Scale Factor/ External Command Sensitivity

Lab Bench     Input Signals       Input Signals     Signals	🖻 Lab Bench 🔁 Scale Factor Assistant: Axopatch 200 series
Analog OUT #0 Analog OUT #1 Analog OUT #2 Analog OUT #3 Cmd 0 L clamp OUT 0 V_clamp VC Cmd Cmd 0 L clamp OUT 0 V_clamp VC Cmd	With some extra information about the state of your amplifier the Scale Factor Assistant will automatically choose the command scale factor. Fill in the form below to have the command scale factor automatically calculated: 1. Mode Setting Specify the setting of the Mode switch on the front panel of the Axopatch: Track C V-Clamp
Signal units: m 💌 V Default	C I=0 C I-Clamp Normal C I-Clamp Fast
Scale factor (mV/V):     100     Scale Factor       Holding level     AultiClamp 700B (Demo)     Image: Comparison of the sector of the s	2. Config Setting     Specify the setting of the Config switch on the front panel of     the Axopatch:     Patch     Patch     Whole Cell (B=1)     Whole Cell (B=0.1)
Quick Select     Advanced     About       General     Gains     Auto     Audio       Image: Channel 1     Channel 2	3. Ext. Command Input Specify which external command input you have the command voltage connected to:
Voltage Clamp Feedback Resistor Experiment Type Range C 50 MΩ Whole Cell 1 - 200 nA C 500 MΩ Whole Cell 0.1 - 20 nA	C 20 mV/V ● 100 mV/V (Axopatch 200B only) Scale factor = 100 mV/V.
C         5 GΩ         Patch         10 - 2000 pA           C         50 GΩ         Patch         0.2 - 200 pA           External Command Sensitivity	Help < Back Finish Cancel
C OFF C 20 mV/V ● 100 mV/V	

Molecular Devices

#### Together through life sciences.

## **Scale factor/Output**

🖻 Lab Bench 🛛 🔀	🔁 Edit Protocol - (untitled)
Input Signals Output Signals Digitizer Channels Signals Analog OUT #0 Analog OUT #1 Analog OUT #2 Analog OUT #3 Cmd 0 Cmd 0 Cmd 0 Cmd 0 L clamp OUT 0 V_clamp VC Cmd Rename VC Cmd	Mode/Rate       Inputs       Outputs       Trigger       Statistics       Comments       Math       Waveform       Stimulus         Analog OUT Channels
Cmd 0       Default         Signal units:       m       ✓       Default         Scale factor (mV/V):       100       Scale Factor         Holding level (mV):       -70       Assistant	

#### Together through life sciences.



# Q & A

Q: I want to perform current clamp experiment. How do I change the command voltage to current


Analog OUT Channels				
Channel #0:	Cmd 0	► Range (mV): -200.00 to 199.99 at 20.0 mV/V		

 A: You need to change the signal unit and scale factor on the Lab Bench.





.ab Bench		
Input Signals Output Si Digitizer Channels	gnals   _ Signals	
Analog OUT #0 Analog OUT #1 Analog OUT #2 Analog OUT #3	Cmd 0 Cmd 0 I clamp IClamp1 OUT 0 V_clamp VC Cmd	Add Delete Rename
- I_clamp Signal units: _ p _	• A	Default
Scale factor (pA/V): Holding level (pA):	400	Scale Factor Assistant



#### Together through life sciences.



# Trigger

P	🖀 Edit Protocol - (untitled)						
	Mode/Rate Inputs Outputs Trigger S						
	Start trial with:						
	Immediate 🗨						
	Tr Immediate Digitizer START Input Space Bar Line Frequency						
	Software Message (Start Only) – Software Message (Start and Stop)						

### Immediate

 As soon as you choose <u>Acquire / Record</u> or <u>Acquire / View Only</u>, or press the equivalent toolbar button, the trial starts.

### Digitizer START Input

- To trigger the start of a trial from an external device.
- Space Bar
  - Start the trial with a toolbar press or mouse click.
- Line Frequency
  - Only for series 132x Axon digitizers
  - synchronized with the mains line frequency (typically 50 or 60 Hz).
- Software Message
  - To trigger the trial from the other program



#### Together through life sciences.

# Trigger

3	Edit Protocol - (untitled)					
	Mode/Rate Inputs Outputs Trigger 9					
	Start trial with:					
	Software Message (Start and Stop) 💌					
	Trigger source:					
	Internal Timer					
	Internal Timer					
	Digitizer START Input					
	Space Bar Line Frequency					
	End frequency					

### Internal Timer

 Determines acquisition according to the length of the sweeps and Start-to-Start Intervals and any Conditioning Trains or P/N Leak Subtraction

### First Acquired Signal

• Once a threshold in the signal from the first Analog IN Channel is crossed.

### Digitizer START Input

- To trigger the start of a trial from an external device.
- Space Bar
  - Acquisition is started by pressing the space bar
- Line Frequency
  - Only for series 132x Axon digitizers
  - Synchronized with the mains line frequency (typically 50 or 60 Hz).



#### Together through life sciences.

# Scope Trigger


🔽 Scope trigger				
Statistics Settings				
Pretrigger length (ms):	5	= 50 samples	Polarity Rising	
Statistics threshold (pA):	4780		C Falling	
Hysteresis				

- A 5 V TTL trigger signal from a digitizer output when the signal reaches the threshold.
  - The rear panel TRIGGER OUTPUT BNC (Digidata 1322)
  - The front panel SCOPE BNC (Digidata 1440)





### **Scope trigger**






Molecular Devices

# Hysteresis


✓ Scope trigger Statistics Settings Hysteresis		
Pretrigger length (m Air A Time hysteresis (ms):	0.1 = 1 samples	
Statistics threshold	0.10 = 19.84 pA	
Hysteresis		

- To prevent signal noise activating false triggers, you can adjust the hysteresis settings to reduce trigger sensitivity.
- Time hysteresis
  - This field adjusts the amount of time that the signal has to be under the threshold level (as determined by the Level Hysteresis setting below) to re-arm the trigger.
- Level hysteresis
- This field adjusts the distance that the signal must be under the threshold in order to re-arm
   Together the threefetrigger (subject to the Time Hysteresis
   Collector research use thing a brook of the median are property of Molecular Devices, LLC or their respective owners.





## **Time hysteresis**



 Adjust the time hysteresis to reduce sensitivity to avoid the false triggers.





#### Together through life sciences.

## **Threshold-Based Statistics**

✓ Threshold-Based Statistics	
Percentage above threshold Mean open (ms)	
Event frequency (Hz) Mean closed (ms)	
Statistics update period (s): 1 = 10000 samples	
Always save statistics at the end of each recording	ing

### Percentage above threshold

 Display the percentage of time that the signal is above the threshold value


### Event frequency (Hz)

- Display the frequency of threshold-crossing
- Mean open (ms)
  - Display the average time that the signal is above threshold.
- Mean closed (ms)
  - Display the average time that the signal is below threshold.



#### Together through life sciences.

### **Threshold-Based Statistics**








### **Statistics**


🖻 Edit Protocol - (untitled)	×	
Mode/Rate       Inputs       Outputs       Trigger       Statistics         ✓       Shape Statistics       ✓       Analog IN Signal       IN 0       ✓         ✓       Positive-going       Selected signals:        Measurements          ✓       Negative-going       Selected signals:       ✓       ✓       Peak amplitude (pA)         ✓       No       ✓       Time of peak (ms)          ✓       Absolute       ✓       Time of antipeak (ms)         Baseline Region       Cursor Region       ✓       Standard deviation (pA)         ✓       Area (pA · ms)		
From (ms):       0       => sample 1       +       Half-width (ms)         To (ms):       11.1       => sample 112       +       Maximum rise slope (pA/ms)         Time of maximum rise slope (ms)       +       Maximum decay slope (pA/ms)		
✓ Search Region       1       ↓       Time of maximum decay slope (ms)         Range:       Cursor Region       ✓       Slope (pA/ms)         From (ms):       22.2       => sample 223       +       Rise slope (pA/ms)         To (ms):       33.3       => sample 334       Decay slope (pA/ms)       from       10 % to       90 %         Locay time (ms)       30.8       +       Decay time (ms)       from       90 % to       10 %		
Smoothing window (samples): 1 + Clear after saving		
OK Cancel Help Acquisition mode: Episodic stimulation Update Preview		












Devices



# Q & A

- How many search regions can be applied during data acquisition?
- How?

Clampex - Scope - Demo Digitizer
e Edit Acquire View Configure Tools Window Help
j
1456 INBA CALE FAR FI
$\sum_{c_1} C_{c_2} C_{c_3} C_{c_4} C_{c_5} C_{c_6} C_{c_7} C_{c_8} C_{A1} C_{A2} C_{A3} C_{A3}$
Episodic md O (mV) Scope - Demo Digitizer
md 1 (mV) 10000 (0
<u>md 2 (nA)</u> 8000-





### Comments


Edit Protocol - Current Clamp 1.pro	×
Mode/Rate       Inputs       Outputs       Trigger       Statistics       Comments         ✓       Comments       ✓       Comments       ✓         ✓       Write each trial       ✓       Prompt each trial       ✓         Information       Numeric identifier #1 (e.g. temperature):       25       ✓         Numeric identifier #2 (e.g. pressure):       50       ✓         Numeric identifier #3 (e.g. concentration):       100       ✓	
File comment: NMDA receptor	
Amplifier mode: Voltage Clamp	

- The comments can be written automatically to every data file
- The comments are stored with data file information in the file header, viewed from <u>File / Properties</u>, and can also be displayed in the <u>Data File Index</u>.





### Comments



Together through life sciences.



### Math


Edit Protocol - Curre     Mode/Rate   Inputs   Outp     Math Signal     Data Source     Analog IN #A: IN 0	uts Trigger Statistics Comments Math Waveform Stimulus	
Equation © General purpose Operator © Addition © Subtraction © Multiplication © Division	C Ratio dyes         Constants         K1:       1         K3:       1         K4:       0         K5:       0         K6:       0	
Yielding:	+K2) <op>(K3*#B+K4) +0)+(1*10Vm+0)</op>	

- Allows arithmetic calculation between two analog input signals
- Calculate on-line ratios, products, sums and differences.



#### Together through life sciences.

### Math



#### Together through life sciences.



### Waveform

#### X dit Protocol - (untitled) Mode/Rate Inputs Outputs Trigger Statistics Comments Math Waveform Stimulus Waveform Analog OUT: Cmd 0 Info - Analog Waveform Digital Outputs 🖲 Epochs 🔘 Stimulus file Use holding Intersweep holding level: 100 ms Epoch Description В Α 112 mV Step Off Туре Sample rate 92 mV Fast Fast First level (mV) 112 72 mV -20 Delta level (mV) 100 First duration (ms) Delta duration (ms) 0 Digital bit pattern (#3-0) Digital bit pattern (#7-4) 0000 Train rate (Hz) Pulse width (ms) 0 mV Number of sweeps = 10 Stimulus File.. Summary Channel #2 Channel #1 Cha Channel #0 Acquisition mode: ΟK Help Update Preview Cancel Episodic stimulation

#### Together through life sciences.




### **Epochs**

dit Protocol - (untitle	d)					
Mode/Rate   Inputs   Out	puts Trigger Statistics					
Waveform Analog OUT: Cmd 0 Info Analog Waveform  Epochs O Stimulus file Intersweep holding level: Use holding						
	A B C					
Type Sample rate	Step Off Off					
First level (mV)	Off					
Delta level (mV)	<b>√ Step</b> Ramp					
First duration (ms)	Pulse train					
Delta duration (ms)						
Digital bit pattern (#3-0)	Triangle train					
Digital bit pattern (#7-4)	Cosine train					
Train rate (Hz) Pulse width (ms)						

#### Together through life sciences.




### **Default waveforms**



#### Together through life sciences.



### **Default waveforms**




Together through life sciences.



lit Protocol - (untitled		inin Í Comu	nts Math	. Wa	eform S			
Mode/Rate   Inputs   Outp Waveform Analog OUT: ✓ Analog Waveform ← Epochs ← Stimu Intersweep holding lev	Cmd 0In	ío	⊡ Digita ⊡ A	l Outputs	n logic for		ins holding	nfo
Epoch Description Type Sample rate First level (mV) Delta level (mV) First duration (ms) Delta duration (ms) Digital bit pattern (#3-0) Digital bit pattern (#7-4) Train rate (Hz) Pulse width (ms)	A B			F	G	H		
Number of sweeps = 10 Stimulus File Summary Channel #0 Channel #1	file not selected.	Alloc annel #3	ated time:		aveforms	∏ Alter     ■	nate Digiti	al Outputs



#### Together through life sciences.



#### Together through life sciences.



Edit Protocol - (untitled)
fode/Rate   Inputs   Outputs   Trigger   Statistics   Comments   Math   Waveform   Stimulus
Waveform Analog OUT: Cmd 0       Info         Analog Waveform       Digital Outputs         Epochs       Stimulus file         Intersweep holding level:       Use holding
Epoch Description A B C D E F G H I J
Type         Sample rate
First level (mV)
Delta level (mV)
First duration (ms)
Delta duration (ms)
Digital bit pattern (#3-0)
Digital bit pattern (#7-4)
Train rate (Hz)
Pulse width (ms)
Number of sweeps = 10 Allocated time:
Stimulus File Stimulus file not selected.
Summary
Channel #0 Channel #1 Channel #2 Channel #3 Alternate Waveforms Alternate Digital Outputs
OK Cancel Help Acquisition mode: Episodic stimulation Update Preview

#### Together through life sciences.



Stimulus File		? 🔀
Look in: 🗀 2010 Biophysic	al Meeting 📃 🗲 🛍	) 💣 🎟 -
🖺 AP wave2.abf		
File name: AP wave2.abf		ОК
Files of type: Axon Binary Fi	ile (*.abf)	Cancel
Signal Specifications Use Analog IN signal:	IN #15 (ADC #15)	Help
Use sweep number (1-1):	All Sweeps 🗨	
Gain factor:	1	
Offset (mV):	0	



#### Together through life sciences.

#### 🚵 Edit Protocol - (untitled)


Waveform Analog OUT: → Analog Waveform → C Epochs ⓒ Stimu Intersweep holding lev	lus file	li se holdin	nfo	]	M A	<b>I Outputs</b> ctive high weep bit	-	 ains holding	Info
Epoch Description Type Sample rate First level (mV) Delta level (mV) First duration (ms) Delta duration (ms) Digital bit pattern (#3-0) Digital bit pattern (#7-4)	A	B				F	G		
		wave2.ab hannel # el #2   C		eeps,	ed time:	,		mata Dia	yital Outputs



©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.



#### **Action Potential Waveform**







### **Digital Outputs**

#### Edit Protocol - (untitled)

Mode/Rate       Inputs       Outputs       Trigger       Statistics       Comments       Math       Waveform       Stimulus         Waveform Analog OUT:       Cmd 0       Info       Info       Info       Info         ✓       Analog Waveform       ✓       Digital Outputs       Info         ✓       Epochs       Stimulus file       ✓       Active high logic for digital trains       Info         Intersweep holding level:       Use holding       ✓       Intersweep bit pattern:       Use holding       ✓										
Epoch Description	A	В		D	E	F	G	Н		J
Туре	Step	Off								
Sample rate	Fast									
First level (mV)	112	0	0	0	0	0	0	0	0	0
Delta level (mV)	-20	0 🖊	0	0	0	0	0	0	0	0
First duration (ms)	100	2	0	0	0	0	0	0	0	0
Delta duration (ms)	0 🕨	0	0	0	0	0	0	0	0	0
Digital bit pattern (#3-0)	1111	0000	0000	0000	0000	0000	0000	0000	0000	0000
Digital bit pattern (#7-4)	0000	0000	0000	0000	0000	0000	0000	0000	0000	0000
Train rate (Hz)	0	0	0	0	0	0	0	0	0	0
Pulse width (ms)	0	0	0	0	0	0	0	0	0	0
Number of sweeps = 10       Allocated time:       106.2 of 200 ms         Stimulus File       Final level -68.00 mV         First duration 100.00 ms (1000 samples)         Summary										

Together through life sciences.





#### **Digital Bit Pattern**


# Digital bit pattern (#3-0): 0101 $\int_{D0\#3 D0\#2} \int_{D0\#1 D0\#0} \int_{D0\#0} \int_{D0\#1 D0\#0} \int_{D0\#0} \int$

1→ High, single TTL (5V) 0→ Low, No single TTL (5V) an asterisk,\*→ A pulse train of TTL (5V)



#### Together through life sciences.

### **Digital Bit Pattern---Single Pulse**




©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.

Together through life sciences.



### **Digital Bit Pattern---Train Pulse**






©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.

Together through life sciences.

### **Digital Bit Pattern---Train Pulse**




#### Together through life sciences.



### **Digital Bit Pattern**


	A	В	С	D	E	F	G	Н		J
Туре	Step	Pulse	Step	Pulse	Step	Pulse	Step	Pulse	Off	Off
Sample rate	Fast	Fast	Fast	Fast	Fast	Fast	Fast	Fast	Fast	Fast
First level (mV)	0	0	0	0	0	0	0	0	0	0
Delta level (mV)	0	0	0	0	0	0	0	0	0	0
First duration (ms)	1000	1000	1000	1000	1000	1000	1000	1000	0	0
Delta duration (ms)	0	0	0	0	0	0	0	0	0	0
Digital bit pattern (#3-0)	0000	000×	0000	000×	0000	000×	0000	000×	0000	0000
Digital bit pattern (#7-4)	0000	0000	0000	0000	0000	0000	0000	0000	0000	0000
Train rate (Hz)	0	10	0	10	0	10	0	10	0	0
Pulse width (ms)	0	1	0	1	0	1	0	1	0	0
						_				



Together through life sciences.




### **Intersweep Holding Level**





#### Together through life sciences.


## Intersweep Holding Level---Use last epoch



#### Together through life sciences.




## Intersweep Holding Level---Use Holding




#### Together through life sciences.



### **Alternative Waveform**











### **Alternative Waveform**










### **Alternative Waveform**




Together through life sciences.



### **Alternative Digital Outputs**






Together through life sciences.



### **Alternative Digital Outputs**








### **Pre-sweep Train**

- A pre-sweep train consists of repeated square waveform pulses
- Conditioning trains
- No acquisition occurs



#### **Pre-sweep Train**




#### Together through life sciences.




#### **P/N Leak Subtraction**

- Leak subtraction corrects for the passive responses by applying test subsweeps before or after the main stimulus waveform
- The technique is called P/N subtraction, in which N subsweeps each 1/Nth of the amplitude of the main stimulus waveform are applied.
- Example: Number of subsweep (N)=4, the stimulus waveform (P)= 80 mV, the pulse of each subsweep=80/4=20 mV



### **P/N Leak Subtraction**



evices


#### Toaether through life sciences.

### User List

- Customizing output features
  - Analog
  - digital
- Overrides the generalized settings made elsewhere in the Protocol Editor.



### **User List---Holding Level**



## You can set the arbitrary delta values for the waveform holding levels.


#### Together through life sciences.





#### **User List---Epoch Duration**

Sweep 1: 22 ms Sweep 2: 44 ms Sweep 3: 100 ms Sweep 4: 75 ms



## You can set the arbitrary delta values for the waveform sweep durations.

#### Together through life sciences.





#### **User List---Change in Parameters**

🔽 User List		
Parameter to change:	Number of pulses in pre-sweep train	-
List of parameter values:	Number of pulses in pre-sweep train Pre-sweep train baseline duration (ms)	~
	Pre-sweep train baseline level	_
Membrane Test Betwe	Pre-sweep train step duration (ms) Pre-sweep train step level	
Settling time (ms): Before	Post-train duration (pre-sweep train) (ms)	24
	Post-train level (pre-sweep train)	
Channel #0 Channel #1	Channel #2 Channel #3	

- Epoch A–J level
- Epoch A–J duration
- Epoch A–J digital pattern
- Epoch A–J train period
- Epoch A–J train pulse width
- Time between sweep starts (s)
- Inactive analog OUT holding level
- Digital intersweep holding level
- Number of P/N subsweeps

- Number of pulses in pre-sweep train
- Pre-sweep train baseline duration (ms)
- Pre-sweep train baseline level
- Pre-sweep train step duration (ms)
- Pre-sweep train step level
- Post-train duration (pre-sweep train) (ms)
- Post-train level (pre-sweep train)




#### Together through life sciences.

### Lists of parameter values

- The list of values for the <u>A Epoch level</u> might be:
  -50, 65, -20, 100, 29
- The list of values for the <u>A Epoch duration</u> might be:
  22, 44, 100, 75

#### **Membrane Test Between Sweeps**



#### = Membrane test

#### Together through life sciences.





#### **pCLAMP**<sup>TM</sup> 10 data acquisition software: Key features review



Together through life sciences.

### **Key features of Clampex 10**

- Online Statistics
- Sequencing keys
- User list
- Membrane test between sweep
- Protocol editor
- LTP assistant
- Stimulus file
- Triggering external devices



#### Together through life sciences.


### **Online Statistics**

- Measure various parameters of evoked events such as peaks, slopes, areas, and rise time
- Available in episodic and oscilloscope modes
- Up to 8 search regions



#### **Outward signal**


#### **Inward signal**




Example: Excitatory postsynaptic potential (EPSP) Example: Excitatory postsynaptic current (EPSC) Population spikes Field EPSP



#### Together through life sciences.









🛜 Clampex - 10d18005.abf - Scope	- Demo Digitizer
File Edit Acquire View Configure Tools	; Window Help
and the second	Window Help
	OK Cancel Help Acquisition mode: Episodic stimulation Update Preview




#### Together through life sciences.

<ul> <li>Clampex - 10d18005.abf - Scope - D</li> <li>File Edit Acquire View Configure Tools V</li> <li>New Protocol</li> <li>Open Protocol</li> <li>Save Protocol</li> <li>Save Protocol As</li> </ul>	
Edit Protocol Waveform Preview Export protocol ( pCLAMP 9 compa Cmd 0 (mV 0 Cmd 1 (nA 0 0 Pause-View	★ Edit Protocol - Post-LTP and Memb.pro         Mode/Rate       Inputs       Outputs       Trigger       Statistics         ✓ Analog IN Signal       IN 0       ✓       Measurements         ✓ Positive-going       Selected signals:       Measurements         ④ Ø Peak amplitude (pA)       Time of peak (ms)         ● Absolute       Time of peak (ms)         Baseline Region       Cursor Region       Mean (pA)         From (ms):       0.75       => sample 16         To (ms):       1.6       => sample 33         ✓ Search Region       1       Ime of maximum rise slope (ms)         +       Maximum decay slope (ms)       Time of maximum decay slope (ms)         +       Baseline (pA)       Slope (pA/ms)         ✓ Search Region       ✓       Baseline (pA)         From (ms):       1.6       Slope (pA/ms)         Time of maximum decay slope (ms)       Time of maximum decay slope (ms)         •       Baseline (pA)       Baseline (pA)
Click "Help" for feature description	From (ms):       3.15       => sample 64       +       Rise time (ms)       } from       10 % to       90 %         To (ms):       9.55       => sample 192       +       Decay slope (pA/ms)       +       Max       Min         Smoothing window (samples):       1       +       Decay time (ms)       } from       90 % to       10 %         OK       Cancel       Help       Acquisition mode:       Update Preview

### **Sequencing Keys**

- Sequencing keys command allows you associate events, or a sequence of event, with a keystroke.
  - Set various digital outputs
  - Change the holding levels
  - Insert a comment tag
  - Start a Membrane Test
  - Load or run a protocol
  - Display a prompt
- Use sequencing keys to link one event to another, and run an experiment in an "automated" fashion.




#### Molecular Devices

### **Sequencing Keys**

- The key strokes and tool button identifiers used for sequencing keys are:
  - <Alt + 0> through <Alt + 9> (i.e. hold down the Alt key and press the relevant numeral)
  - <Alt + Shift + 0> through <Alt + Shift + 9>
  - <Ctrl + 0> through <Ctrl + 9>
  - <Ctrl + Alt + 0> through <Ctrl + Alt + 9>
  - <Ctrl + Shift + 0> through <Ctrl + Shift + 9>
- 50 sequencing keys can be triggered in one time.





#### Together through life sciences.




### **Sequencing Keys**


 Only the key for the first element needs to be triggered to start off the entire sequence.








#### Together through life sciences.



#### Together through life sciences.









### **User List**

- Provide a way of customizing one of a range of analog and digital output features.
- It overrides the generalized settings made elsewhere in the Protocol Editor.



### User List

- With User List you can select one of a range of parameters offered, to configure on a sweep-bysweep basis.
- This function allows you to set specific values for each sweep in a run. You can also set a sequence of values and have this repeat.



### **Change in Epoch holding levels**

Sweep 1: -50 mV Sweep 2: 65 mV Sweep 3: -20 mV Sweep 4: 100 mV Sweep 5: 29 mV


# You can set the arbitrary delta values for the waveform holding levels.

#### Together through life sciences.



### **Change in Epoch holding levels**

Sweep 1: -20 mV Sweep 2: 100 mV Repeat...








### **Change in Epoch durations**

Sweep 1: 22 ms Sweep 2: 44 ms Sweep 3: 100 ms Sweep 4: 75 ms



# You can set the arbitrary delta values for the waveform sweep durations.

#### Together through life sciences.



### **Change in Epoch durations**

Sweep 1: 22 ms Sweep 2: 75 ms Repeat...



Molecular Devices

### **Parameter to Change**

- Epoch A–J level
- Epoch A–J duration
- Epoch A–J digital pattern
- Epoch A–J train period
- Epoch A–J train pulse width
- Time between sweep starts (s)
- Inactive analog OUT holding level
- Digital intersweep holding level
- Number of P/N subsweeps

- Number of pulses in pre-sweep train
- Pre-sweep train baseline duration (ms)
- Pre-sweep train baseline level
- Pre-sweep train step duration (ms)
- Pre-sweep train step level
- Post-train duration (pre-sweep train) (ms)
- Post-train level (pre-sweep train)



#### Together through life sciences.

List of parameter values: Number of pulses in pre-sweep train Pre-sweep train baseline duration (ms) Pre-sweep train step duration (ms) Pre-sweep train step duration (ms) Pre-sweep train step level Settling time (ms): Before Post-train duration (pre-sweep train) (ms) Post-train level (pre-sweep train)	ode/Rate   Inputs   Out	tputs Trigger Statistics Comments Math Waveform Stimulus	
Number of pulses in train:  Pulse frequency = 500 Hz Train duration = n/a Baseline duration (ms):  Step duration (ms):  Post-train duration (ms):  Post-train level (mV):  Post-train level (mS) Post-train level (mV):  Post-train level (mV): Post-train level (mV): Post-train level (mV): Post-tra		: Cmd 0 Info	
Baseline level (mV):       0       Step level (mV):       0       Post-train level (mV):       0         P/N Leak Subtraction       Apply to Analog IN signal:       IN 0       Execution       Polarity       Same as waveform         Number of subsweeps:       4       Image: Composite to waveform       Opposite to waveform       Opposite to waveform         Subsweep start-to-start (ms):       Minimum       Subsweep holding level (mV):       0         Settling time (ms):       100       Image: Show corrected sweep data         Ist of parameter values:       Number of pulses in pre-sweep train         Pre-sweep train baseline duration (ms)       Pre-sweep train baseline level         Pre-sweep train baseline level       Pre-sweep train step level         Setting time (ms):       Before       Image: Pre-sweep train step level         Pre-sweep train baseline level       Pre-sweep train step level       Pre-sweep train step level         Setting time (ms):       Before       Post-train duration (pre-sweep train)       Image: Pre-sweep train level (pre-sweep train)		n: 1 🚔 Pulse frequency = 500 Hz Train duration = n/a	
<ul> <li>P/N Leak Subtraction         Apply to Analog IN signal:         IN 0         IN 0</li></ul>	Baseline duration (ms):	1 Step duration (ms): 1 Post-train duration (ms): 10	
Apply to Analog IN signal: IN 0 Number of subsweeps: 4 Subsweep start-to-start (ms): Minimum Subsweep holding level (mV): 0 Settling time (ms): 100 Settling time (ms): 100 Show corrected sweep data User List Parameter to change: Number of pulses in pre-sweep train Parameter values: Number of pulses in pre-sweep train Pre-sweep train baseline duration (ms) Pre-sweep train baseline level Pre-sweep train step level Settling time (ms): Before Polarity Repeat Pre-sweep train step level Settling time (ms): Before Post-train duration (pre-sweep train) (ms) Pre-sweep train level (pre-sweep train) (ms)	Baseline level (mV):	0 Step level (mV): 0 Post-train level (mV): 0	
Apply to Analog IN signal: IN 0 Setting time (ms): Minimum Subsweep holding level (mV): 0 Setting time (ms): 100 Show corrected sweep data User List Parameter to change: Number of pulses in pre-sweep train List of parameter values: Number of pulses in pre-sweep train Pre-sweep train baseline duration (ms) Pre-sweep train step duration (ms) Pre-sweep train step level Setting time (ms): Before Post-train duration (pre-sweep train) Post-train level (pre-sweep train) Post-train level (pre-sweep train)	P/N Leak Subtraction		
Subsweep start-to-start (ms): Settling time (ms): User List Parameter to change: List of parameter values: Mumber of pulses in pre-sweep train Pre-sweep train baseline duration (ms) Pre-sweep train baseline level Pre-sweep train step duration (ms) Pre-sweep train step level Settling time (ms): Before Settling time (ms): Before Post-train duration (pre-sweep train) Pre-sweep train (ms) Pre-sweep train step level Post-train duration (pre-sweep train) Pre-sweep train (ms) Pre-sweep train (ms) Pre-sweep train step level Post-train duration (pre-sweep train) Pre-sweep train (ms) Pre-sweep trai	Apply to Analog IN sign	-6 010	
Setting time (ms):       100       Image: Show corrected sweep data         Image: User List       Parameter to change:       Number of pulses in pre-sweep train         List of parameter values:       Number of pulses in pre-sweep train         Image: Pre-sweep train baseline duration (ms)       Image: Pre-sweep train step duration (ms)         Pre-sweep train step duration (ms)       Pre-sweep train step level         Setting time (ms):       Before         Post-train duration (pre-sweep train)       Image: Post-train level (pre-sweep train)	Number of subsweeps:	4 C After C Opposite to waveform	
<ul> <li>User List</li> <li>Parameter to change:</li> <li>List of parameter values:</li> <li>Mumber of pulses in pre-sweep train</li> <li>Pre-sweep train baseline duration (ms)</li> <li>Pre-sweep train step duration (ms)</li> <li>Pre-sweep train step level</li> <li>Settling time (ms): Before</li> </ul>	Subsweep start-to-start	(ms): Minimum 💽 Subsweep holding level (mV): 0	
<ul> <li>✓ User List</li> <li>Parameter to change:</li> <li>List of parameter values:</li> <li>Number of pulses in pre-sweep train Pre-sweep train baseline duration (ms) Pre-sweep train baseline level</li> <li>Membrane Test Betwy Pre-sweep train step duration (ms)</li> <li>Pre-sweep train step level</li> <li>Post-train duration (pre-sweep train)</li> </ul>	Settling time (ms):	100 Show corrected sweep data	
1 Anna 1 40   Elemented 11   Elemented 12   Elemented 12	Viser List Parameter to change:	Number of pulses in pre-sweep train	
	Membrane Test Bet Settling time (ms): Befor	<ul> <li>Pre-sweep train baseline duration (ms)</li> <li>Pre-sweep train baseline level</li> <li>Pre-sweep train step duration (ms)</li> <li>Pre-sweep train step level</li> <li>rane Test Set</li> </ul>	

#### Together through life sciences.




### **Examples: List of Parameter values**

- The list of values for the <u>A Epoch level</u> might be:
  -50, 65, -20, 100, 29
- The list of values for the <u>A Epoch duration</u> might be:
  22, 44, 100, 75




### **Change in Epoch holding levels**



Together through life sciences.



### **Change in Epoch durations**



Together through life sciences.



### **Membrane Test Between Sweeps**

Together through life sciences.



### **Membrane Test**


- Access resistance, Ra
- Membrane resistance, Rm
- Membrane capacitance, Cm

- Time constant, Tau
- Holding current, Hold




Edit Protocol - Post-LTP protocol.pro		×
Mode/Rate Inputs Outputs Trigger Statistics Comments Ma	ath   Waveform   Stimulus	1
Acquisition Mode         File       C Gap-free         C       Gap-free         C       Variable-length events         C       High-speed oscilloscop	<ul> <li>Episodic stimulation</li> </ul>	
Trial Hierarchy       Trial delay (s):       0       Runs/trial:       1       Sweeps/run:       100	Sweep (s): 1	h Waveforr Stimulus
Epis         Sweeps/run.         100         = 3.01 MB           Cmd I         Sweep duration (s):         1         = 20000 samples           0         First noiding:         E pochs:         Last holding:           Cmd I         15.6 ms         968.8 ms         15.6 ms           0         312 samples         19376 samples         312 samples		ain duration = 2 ms Post-train duration (ms); 10
Cmd 1         Sampling Rate per Signal           0         Fast rate (Hz): 20000           Cmd ( 0         Fast rate (Hz): 20000           0         Slow rate (Hz): 20000	Averaging Runs/trial = 1, no averaging Options	Post-train level (mV): 0 Polarity Same as waveform
Space available is 24114 sweeps = 50071 MB Total data Allow automatic analysis in other programs	throughput is 20 kHz (= 2.29 MB/min)	C Opposite to waveform ng level (mV): ted sweep data
OK Cancel Help Acquisition Episodic	c stimulation	
	Test Between Sweeps s): Before 100 After 100 For oth	er settings, use Configure / Membrane Test Setup
Channel #0 Ch OK	nannel #1 Channel #2 Channel #3 Cancel Help Acquisit	ion mode: odic stimulation Update Preview
gether through life sciences. 012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of M		Molecula Device

#### Together through life sciences.



#### Together through life sciences.





©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.

#### Devices

🖞 Cla	mofit	- [Resi	ılts1 - F	Results]													
_		-			Event Detect	tion Confi	igure Tools	Window	Help								
	New New Dat	a File In	ndex														
	Open Da Open Da		ons			Ctrl+0											
	Open Da	•					BR. A	Xal	III Kas								
É	Open Ot	her					Lab Be		Ł								
	Close					Ctrl+W	Resul	:s & Statist		l <mark>ampfit - [Res</mark> ile Edit View			tection Configu	ure Tools Window	Help		
	Close All								<u> </u>			ì 🕹 🖻					
									Factor References		「「「「「「」」						
									1	<b> </b> <u>∧</u> 4 X   Y+	XY t	<u> </u> <u> </u> <u> </u> <u></u>	≣ <u>≫</u> Σ	Mu Ac ≚c #	ų <b>t</b>		
														Memb Test 0-Mem		Memb Test M	lemb Test
									1								
									2		-396.618	-341.907	0.0244	1023.98	640.185	0	15.625
									4	1.00001	-392.27	-304.071	0.0244	1023.98		0	15.625
									5	2.00001	-425.617	-327.669	0.0244	1023.98		0	15.625
									6	3.00001	-389.135	-369.928	0.0244	1023.98	640.185	0	15.625
									7	4.00001	-401.418	-527.35	0.0244	1023.98	640.185	0	15.625
									8	5.00001	-407.63	-301.137	0.0244	1023.98		0	15.625
									9 10 11	6.00001	-377.328	-433.741	0.0244	1023.98		0	15.625
									10	7.00001	-429.17	-300.389	0.0244	1023.98		0	15.625
									11	9.00001	-380.015 -407.305	-241.361 -269.556	0.0244	1023.98		0	15.625
									12	9.00001	-407.305	-269.556	0.0244	1023.98		0	15.625
									13		-410.625	-317.562		1023.98		0	15.625
									14			-402.452	0.0244	1023.98		0	15.625
									16		-411.306	-402.432	0.0244	1023.98		0	15.625
									17		-406.689	-337.691	0.0244	1023.98		0	15.625
									18			-291.797	0.0244	1023.98		0	15.625
									19		-411.435	-325.913		1023.98		0	15.625
									20		-382.614	-429.453	0.0244	1023.98		0	15.625
									21		-396.152	-312.506	0.0244	1023.98		0	15.625
									22		-392.563	-363.627	0.0244	1023.98	640.185	0	15.625

36

37 38

39

40

41

42

-387.604

-379.512

-416.414

-412.133

-407.789

-387.002

-390.243

-406.754

-401.839

-395.766

-412.66

-372.564

-399.646

-397.646

-401.371

-418.472

-384.828

-377.506

-427.499

-406.899

-260.169

-369.29

-416.725

-336.846

-392.96

-255.915

-214.245

-407.541

-342.663

-317.118

-351.936

-304.056

-371.95

-307.246

-330.406

-330.846

-367.603

-380.116

-278.578

-393.753

0.0244

0.0244

0.0244

0.0244

0.0244

0.0244

0.0244

0.0244

0.0244

0.0244

0.0244

0.0244

0.0244

0.0244

0.0244

0.0244

0.0244

0.0238

0.0238

0.0238

20

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35

36

37

38

39

1023.98

1023.98

1023.98

1023.98

1023.98

1023.98

1023.98

1023.98

1023.98

1023.98

1023.98

1023.98

1023.98

1023.98

1023.98

1023.98

1023.98

995.964

995.964

995.964

640.185

640.185

640.185

640.185

640.185

640.185

640.185

640.185

640.185

640.185

640.185

640.185

640.185

640.185

640.185

640.185

640.185

654.885

654.885

654.885

0

0

0

0

0

0

0

0

0

0

0

0

0

0

0

0

0

0

0

 $\sim$ 

15.625

15.625

15.625

15.625

15.625

15.625

15.625

15.625

15.625

15.625

15.625

15.625

15.625

15.625

15.625

15.625

15.625


15.625

υevices

15.625


аг

Together through life sciences.



Devices

### Summary



©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.

Devices


### **Stimulus File**

 Stimulus Files are used when the Waveform Editor cannot create the desired waveform.

#### Together through life sciences.




#### Edit Protocol - (untitled)



х

evices


### **Default waveforms**



Together through life sciences.




### **Default waveforms**



#### Together through life sciences.



# How to create an action potential waveform?



#### Together through life sciences.






## Advanced Tips for Protocol Writing in the Clampex Data Acquisition Module




#### Together through life sciences.

### Whole cell and Extracellular Recordings





#### Together through life sciences.




#### **Application: Synaptic plasticity, release probability in brain slice preparation**



#### Together through life sciences.

### **Pair-pulse Stimulation: Presynaptic afferents**






#### Together through life sciences.

### **Pair-pulse stimulation**

🔁 Edit Protocol - (untitle	:d)									į	4	🔁 Wa	veform	n Preview				×
Mode/Rate   Inputs   Out	puts   T	rigger   S	tatistics	Comme	nts Ma	th Wa	/eform	Stimulus						1	₩2	3	₹4	
Waveform Analog OUT Analog Waveform Epochs C Stirr Intersweep holding le	Cmd 0 ulus file		Info		Digit	al Output: Active hig sweep bi	s	or digital		Info		Cmd 1 Cmd 0 (mV) (mV)	200 . ▼ 0 - 200 . 200 .	Wave Start				
Epoch Description	A	В	С	D	E	F	G	Н		J		55						
Туре	Step	Step	Step	Step	Step	Off	Off	Off	Off	Off			-200 -					
Sample rate	Fast	Fast	Fast	Fast	Fast	Fast	Fast	Fast	Fast	Fast			200 .					
First level (mV)	-80	-80	-80	-80	-80	0	0	0	0	0		23	0 -					
Delta level (mV)	0	0	0	0	0	0	0	0	0	0		Cmd 2 (m V	- 0				Wave E	nd
First duration (ms)	100	1	49	1	50	0	0	0	0	0			-200					
Delta duration (ms)	0	0	0	0	0	0	0	0	0	0			200					
Digital bit pattern (#3-0)		0001	0000	0001	0000	0000	0000	0000	0000	0000			-					
Digital bit pattern (#7-4)	0000	0000	0000	0008	0000	0000	9999	0000	0000	0000		S da G md S	0-		50 ms			
Train rate (Hz)	0	0	0	0	0	0	0	0	Ũ	0		00						
Pulse width (ms)	0	0	0	0	0	0	0	0	0	0			-200 U					
Number of sweeps = 10				Alloc	ated time		21	0.2 of 3	00 ms				- 0					#
Stimulus File First du	ration 49	).00 ms (4	90 sampl	es)									0					╡
Sumaias Filo																		#
Summary												<u>lĕ</u> ≩						
· · · ·													(	)	100	200	'	300
Channel #0 Channel #1	Chan	nel #2 (	Channel ‡	‡3		ternate W	'aveform	s 🗖 Ali	ternate Dig	gital Outputs				4	Time (m	s) :	Sweep:1 Visible	e:1 of 1 ▶
OK Cano	el	He	elp	Ac	quisition r Episodic	node: : stimulatio	n		Updat	e Preview				=> Clipboard	Close	Help		



#### Together through life sciences.



#### **Application: Synaptic plasticity, LTD in brain slice preparation**



### **Train stimuli-10 pulses**

		🔁 Waveform Preview	×
🚵 Edit Protocol - (untitled)	<u>×</u>		
Mode/Rate Inputs Outputs Trigger Statistics Comme	ents Math Waveform Stimulus	200	<u> </u>
Waveform Analog OUT: Cmd 0 Info ✓ Analog Waveform ← Epochs ← Stimulus file Intersweep holding levet: Use holding ▼	<ul> <li>Digital Outputs</li> <li>Active high logic for digital trains Info</li> <li>Intersweep bit pattern: Use holding </li> </ul>	P     0       -200     Wave Start       200     0	
		55	
Epoch Description A B C D	E F G H I J	-200	
Type Step Pulse Step Off	Off Off Off Off Off Off	200 .	
Sample rate Fast Fast Fast Fast	Fast Fast Fast Fast Fast		
First level (mV) -80 -80 0	0 0 0 0 0 0		
Delta level (mV) 0 0 0	0 0 0 0 0 0	UVave F	End
First duration (ms) 100 1000 0	0 0 0 0 0 0	-200 j	
Delta duration (ms) 0 0 0	0 0 0 0 0 0	200 .	
Digital bit pattern (#3-0) 0000 000* 0000	0000 0000 0000 0000 0000 0000		
Digital bit pattern (#7-4) 0000 0000 0000 0000	0000 0000 0000 0000 0000 0000		
Train rate (Hz) 0 10 0 0	0 0 0 0 0	<b>10 Pulses</b>	
Pulse width (ms) 0 1 0 0	0 0 0 0 0 0		
Number of sweeps = 10 Alloc	ated time: 1246.8 of 1500 ms		≠ ∥
Stimulus File First duration 1000.00 ms (10000 samples)	Pulse count 10 Train rate 10.00 Hz (1000 samples) Pulse width 1.00 ms => 1 % (10 samples)	>     > <td></td>	
Channel #0 Channel #1 Channel #2 Channel #3	Alternate Waveforms 🗍 Alternate Digital Outputs	Time (s) Sweep:1 Visib	le:1 of 1
OK Cancel Help Ac	quisition mode: Episodic stimulation	=> Clipboard Close Help	1.



#### Together through life sciences.

### Train stimuli-50 pulses

🖀 Edit Protocol - (untitled)	X	🖀 Waveform Preview							
Mode/Rate Inputs Outputs Trigger Statistics Comm	ents Math Waveform Stimulus		<b>A</b>						
Waveform Analog OUT: Cmd 0 Info									
Analog Waveform	🔽 Digital Outputs	557							
€ Epochs ○ Stimulus file	Active high logic for digital trains	-200 Wave Start	_						
		200							
Intersweep holding level: Use holding	Intersweep bit pattern: Use holding 💌		-						
A B C D	E F G H I J								
Type Step Pulse Step Off	Off Off Off Off Off	-200 1 200 1	-						
Sample rate Fast Fast Fast Fast	Fast Fast Fast Fast Fast								
First level (mV) -80 -80 0	0 0 0 0 0	B   0     Wave End	- 1						
Delta level (mV) 0 0 0	0 0 0 0 0	US Wave End							
First duration (ms) 100 1000 0	0 0 0 0 0	-200 1							
Delta duration (ms) 0 0 0	0 0 0 0 0	200 .	-						
Digital bit pattern (#3-0) 0000 000* 0000 0000	0000 0000 0000 0000 0000 0000								
Digital bit pattern (#7-4) 0000 0000 0000 0000	0000 0669 0000 0000 0000 0000		-						
Train rate (Hz) 0 50 0 0		50 Pulses							
Pulse width (ms) 0 1 0 0	0 0 0 0 0		_						
Number of sweeps = 10 Allo	cated time: 1246.8 of 1500 ms		-						
Stimulus File   First duration 1000.00 ms (10000 samples	s) Pulse count 50		-						
oundros mess	Train rate 50.00 Hz (200 samples)								
Summary	Pulse width 1.00 ms => 5 % (10 samples)		-						
		■■ 0 0.4 0.8 1.2							
Channel #0 Channel #1 Channel #2 Channel #3	🔲 Alternate Waveforms 🔲 Alternate Digital Outputs	Time (s) Sweep:1 Visible:1	of 1						
OK Cancel Help Ac	cquisition mode: Episodic stimulation	=> Clipboard Close Help	//.						






### Summary

- Telegraphing/Lab Bench settings
- Writing protocols for
  - Pair-pulse stimulation/Digital Output
  - Train stimulation



### **Tetanus Stimulation protocol**




#### **Application: Synaptic plasticity, LTP/LTD in brain slice preparation**



#### Together through life sciences.

### Tetanus Stimulation and Postsynaptic Depolarization



Molecular Devices


Together through life sciences.

### **Tetanus Stimulation and Postsynaptic Depolarization**

🔁 Edit Protocol - (untitleo	l)									[	×	🔁 Wa	veform	n Preview						×
Mode/Rate   Inputs   Outp	uts Í Tri	ager Í S	tatistics [	Comme	nts Í Mat	h Wa	veform	Stimulus						1	$\overline{}$	2	3		4	
Waveform Analog OUT: ☐ Analog Waveform — ⓒ Epochs ⓒ Stimu Intersweep holding lev	Cmd 0 Ilus file		Info	,	Digit.	al Output Active hig		or digital I		Info		Cmd 1 Cmd 0 (mV) (mV)	200 . 0 - -200 - 200 .	Wave Star	「			Π		
Epoch Description	A	B		L D	F	F	G	Гн				[5 5	-							
Туре	Step	Step	Step	Step	Step	Step	Step	Step	Step	Off			-200							
Sample rate	Fast	Fast	Fast	Fast	Fast	Fast	Fast	Fast	Fast	Fast			200 .							
First level (mV)	-80	0	-80	þ	-80	0	-80	0	-80	0		22	_							
Delta level (mV)	0	0	0	þ	0	0	0	0	0	0		Cmd 2 (m d 2	0-						N	ave End
First duration (ms)	100	1000	10000	000	10000	1000	10000	1000	100	0			-200							
Delta duration (ms)	0	0	0	þ	0	0	0	0	0	0		<u> </u>	200							
Digital bit pattern (#3-0)	0000	000×	0000	)00×	0000	000×	0000	000×	0000	0000										
Digital bit pattern (#7-4)	0000	0000	0000	0000	0000	0000	0000	0000	0000	0000		Cmd3 (md3	0-							
Train rate (Hz)	0	100	0	100	0	100	0	100	0	0		55	1							
Pulse width (ms)	0	1	0		0	1	0	1	0	0			-200							
Number of sweeps = 10				Alloc	ated time:		354	50 of 40	1000 ms				U U							
Summary			1000 sam									Appebbbbbbbbbbbbbbbbbbbbbbbbbbbbbbbb			,	1	20	,	,	40
Channel #0 Channel #1	Chann	el #2   (	Channel #	3]		ternate W	/aveform:	s 🗌 Alt	ernate Di	igital Outputs				<u>र</u>		T	ime (s)		Sweep:1	Visible:1 of 1
OK Cance	9	He	elp	Ac	quisition n Episodic		on		Updal	te Preview				=> Clipb	oard	Close		Help		



#### Together through life sciences.




## **Application: Synaptic plasticity, LTP/LTD in brain slice preparation and culture cell**



#### Together through life sciences.


### **Dual Whole-cell Patch-clamp Recordings**





#### Together through life sciences.


### **Dual Whole-cell Patch-clamp Recordings**








### Spike timing protocol: EPSP precedes Action Potential





#### Together through life sciences.


### Spike timing protocol: EPSP precedes Action Potential



#### Together through life sciences.




### Spike timing protocol: EPSP precedes Action Potential



#### Together through life sciences.



## Spike timing protocol: Action Potential precedes EPSP






©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.

Together through life sciences.


# Spike timing protocol: Action Potential precedes EPSP



#### Together through life sciences.



## Spike timing protocol: Action Potential precedes EPSP



#### Together through life sciences.





#### **Action Potential Analysis in Clampfit**



#### Together through life sciences.

#### **Frequent Ask Questions in Action Potential Analysis**

- How can I measure the frequency of action potentials?
- The baseline of my recording is drifting. It is hard to analyze the spikes.
- How do I sort the spikes?
- My recording is too noisy. It is hard to retrieve the spikes from the noises?
- When I do the event search, how can I avoid the large spike of stimulation artifacts?
- How do I perform a phase plot of dV/dt vs V in Clampfit?
- How do I calculate the cardiac action potential duration at 90% repolarization (APD<sub>90%</sub>)?



#### Together through life sciences.

### What we learn today

- Baseline adjustment
- Event Search
- Event sorting
- Noise/Event rejection
- Spike alignment
- Combine trace
- Phase plot
- Action potential analysis



#### Together through life sciences.

### **Features in Clampfit**


- Baseline adjustment
  - Manual baseline adjustment
- Search event
  - Event Detection/Threshold Search
- Event sorting
  - Event Detection/Threshold Search
- Spike alignment
  - Time shift
- Phase plot/Combine trace
  - Arithmetic
- Action potential analysis
  - Statistics

Together through life sciences.






### **Baseline adjustment**







### **Event Search**




If an event crosses the trigger level, the event will be accepted.

#### Together through life sciences.



#### **Event Sorting**

 Extract events from a gap-free file and transform them into sweeps.






#### Together through life sciences.

### **Noise/Event Rejection**

- Event rejection based on amplitude
- Noise rejection based on short duration
- Event rejection based on event length




### **Noise/Event Rejection based on amplitude**




If an event crosses the rejection level, the event will be rejected.

#### Together through life sciences.





### **Noise/Event Rejection based on too short duration**




Noise duration in ms can be set to reject noise spikes.

#### Together through life sciences.



### **Noise/Event Rejection based on event duration**




You can select Min allowed duration to eliminate short events. You can select Max allowed duration to discard events that are too long.

#### Together through life sciences.




### **Peak Alignment**





Molecular Devices

### **Combine traces**



Trace in File A

Together through life sciences.


Trace in File B

#### Trace combined



#### **Phase Plot**

• Plot dV/dt versus V






#### Together through life sciences.


### **Cardiac Action Potential Analysis**

- Resting membrane potential (RMP)
- Action potential Amplitude (APA)
- Action Potential Duration (APD)
- Max Slope of depolarization (Vmax)





### **Event Detection/Threshold Search**





#### Together through life sciences.

#### Set category

- Each category represents a unique search criterion. There are up to 8 different categories
- The positive- and negativegoing describe the location of peak relative to baseline.
  - Action potential is a positivegoing event


Threshold Search	×
Category 1 🐥 Reset /	
<u>Positive-going</u> <u>Negative-going</u>	



#### Together through life sciences.

### **Bring markers**


• Bring markers








#### **Re-arm**




### **Event Rejection**




#### Together through life sciences.




### **Noise Rejection**






Event Start and End without Noise rejection

Together through life sciences.





Event Start and End with 5 ms Noise rejection

#### Molecular Devices

#### ©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.


6/20/2013 Pg. 179

### **Pre-trigger and Post-trigger**




Pretrigger and post-trigger lengths						
🔽 Pretrigger (ms):	5					
🔽 Post-trigger (ms):	20					

Molecular Devices




#### Together through life sciences.

### **Event Rejection**



#### Together through life sciences.




#### Together through life sciences.





### **Event States**



	Trace	Search	Category	State
1	1	1	1	A
2	1	1	1	A
3	1	1	1	AT
4	1	1	1	AT
5	1	1	1	AT


	Trace	Search	Category	State
1	1	1	1	A
2	1	1	1	A
3	1	1	1	S
4	1	1	1	S
5	1	1	1	AT
6	1	1	1	AT
7	1	1	1	AT

Event Detection/Show Suppressed Events



#### Together through life sciences.

### **Phase Plot**



#### Together through life sciences.



### **Summary**

- Baseline adjustment
- Event Search
- Event sorting
- Noise/Event rejection
- Spike alignment
- Combine trace
- Phase plot
- Action potential analysis



#### Together through life sciences.











#### Analysis of Synaptic Events with the Clampfit Data Analysis Module




#### Together through life sciences.

### Agenda

• Postsynaptic miniature currents (EPSCs)

• Pre-drug



• Post-drug

- Are there differences in the amplitude distribution?
  - Post-synaptic mechanism
- Are there differences in the frequency?
  - Pre-synaptic mechanism

#### Together through life sciences.

©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.




abf


## Strategy

# Create a template from one of the data files

Perform a template search on a data segment before addition of compound...




Together through life sciences.






## Strategy

# Create cumulative histograms of the peak amplitudes...



#### ...and the inter-event intervals.







#### Together through life sciences.

## Strategy

Perform a statistical analysis to determine whether there are differences between the two datasets.







#### Together through life sciences.


## **Clampfit Features Discussed**

- Event Detection > Create Template
- Event Detection > Template Search
- Move Cursors efficiently
- View > Zoom > Between Cursors
- Analyze > Arithmetic
- Format > Column > Rename
- Format > Rename Sheet
- Analyze > Histogram
- Efficiently creating a graph using X-Y pairs
- Editing a graph
- Analyze > Kolmogorov-Smirnov Test
- Analyze > Basic Statistics



#### Together through life sciences.


## **Creating a Template**






Together through life sciences.

### **Creating a Template**




#### Together through life sciences.

## **Delimit the Analysis Region**




Full-scale the time axis,



unlock cursors 1 and 2,



#### and send cursor 1...



#### ..."**Home**" to the start of the trace



#### Together through life sciences.


## **Delimit the Analysis Region**

	Double- click on Cursor 2,
	Cursor Properties for Cursor 2         Vertical Cursor Time Value         Move To            • Time (s):         • Sample number:         • 960404         • 10
cond it to	O Tag number:       1       ↓
send it to 30 s,	Number of decimal places to show:       Default (4)         Auto Scale       -         Full Scale       -         Zoom       -         Invert       -
	and zoom the region we are doing to analyze

#### Together through life sciences.




### **Template Search**






#### Together through life sciences.

### **Template Search**



# Define the measurement region.




# Set taus and analysis region.



#### Together through life sciences.

### **Template Search**



#### Together through life sciences.


## **Defining "Sanity-check" Graphs**

✓ ✓ €	Define Graphs
	✓ Graph 1         ○ Conventional histogram       Bin width (ms):       0.1       ✓ Auto         ○ Logarithmic histogram       Bins/decade:       10       ✓ SqRt N         ○ Scatter plot       Convert bin counts to frequency       ✓ Peak ampli       This graph Indicates the absence of an overall trend (run-up or run-down of
	Image: Conventional histogram       Bin width (pA):       0.1       Image: Auto       Category       Image: Auto       Image
	Graph 3 Conventional histo Conventional hist
	Graph 4 distinct amplitude. Conventional histor Logarithmic histogram Scatter plot OK Car OK Car OK Car
	Peak Amplitude (pA) Time Of Peak (ms)

#### Together through life sciences.




### **Copying the Pre-compound Results to a New Sheet**





#### Together through life sciences.


### **Computing the Absolute Peak Amplitude**





#### Together through life sciences.


### **Creating a Cumulative Peak Histogram**





#### Together through life sciences.

### **Repeat the Template Search After Compound Addition**




new file erases all results on the Events sheet (but not Sheet 19)





### **Post-compound Results**




The post-compound data file after completed analysis with

 Together through life sciences.
 "sanity check graphs".

 ©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.




### **Repeat Post-processing Steps**





#### Together through life sciences.

### **Normalizing the Histogram Values**



Normalize the Bin Count columns on the Histogram sheet

#### Rename the two new columns


. »/ I	<b>Ξ</b>   յԱլ   ≙c	≚c	<b>4</b>		
Histo(cBH)	Norm Pre P	F	G	Н	
1	0	0.00117509			
7	0	0.00822562			
15	Column Ren	ame			
50					
6	Name of colu	mn F:		Previous	
86	Norm Post Pe	eak			
10:	,			Next	
111					_
12: 138	OK	Car	icel	Help	
150	L				-
154	0.222222	0.100304			



©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.

Together through life sciences.

### **Creating a Combined, Normalized Cumulative Histogram**



125

Make the first "Bin Center" column an X column

🗎 Res	ults1 - Resu	lts						
	<u>~</u> x   y	XY	1 2 1 C	<b>≞ </b> ≫/ 5	2   J	հղ [≜c	≚c	
	Bin Center		Bin Center	Histo(cBH)	Norr	n Pre P	Norm Post	
1	9.5	Column(s) 1	8.5	1		0	0.00117509	1
2	10.5	1	9.5	7		0	0.00900560	1
3	11.5	6	10.5	15	C	Res	ults1 - Resu	lts
4	12.5	16	11.5	30	C			
5	13.5	26	12.5	50	0.		$\mathbf{N}_{\mathbf{k}} \mathbf{X} \mathbf{V}_{\mathbf{k}}$	X

Make the "Norm Pre Peak" column a Y column


	Bin Center	His X-Y Pair Bin (				
1	9.5					
2	10.5	1				
3	11.5	6				

Define the two columns as an X-Y pair



#### Together through life sciences.

### **Creating a Combined, Normalized Cumulative Histogram**






#### Together through life sciences.


### **Creating a Combined, Normalized Cumulative Histogram**


Graph Window Properties	×
General       X Axis       Y Axis       Plots       Colors/Fonts         Items to Display       V Axis       X major gridlines         X Axis       Y Y axis       X major gridlines         X Name and units       Name and units       X major gridlines         Tick numbers       Rotated text       Y major gridlines         Frame box       Tick numbers       Y major gridlines	
Titles       Image: My EPSPs         Image: Minor title:       Peak Histogram Before and After Compound         Image: Legend:       Legend         Image: Minor title:       Peak Histogram Before and After Compound         Image: Minor title:       Peak Histogram Before and After Compound         Image: Minor title:       Peak Histogram Before and After Compound         Image: Minor title:       Peak Histogram Before and After Compound	

Graph Window Properties	×
General XAxis YAxis Plots Colors/Fonts	
Axis title: Peak Amplitude Units: pA	
Graph Window Properties	X
General X Axis Y Axis Plots Colors/Fonts	
Axis title: Cumulative Fraction	

#### Modify Graph titles, axis titles and plots for an appearance as shown











#### Together through life sciences.

### **Kolmogorov-Smirnov Test**



Molecular Devices


### **Kolmogorov-Smirnov Test Results**







### **Creating Histograms for Inter-event Intervals**






©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.

Together through life sciences.

### **Creating Histograms for Inter-event Intervals**





©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners

Together through life sciences.


### Kolmogorov-Smirnov Test for Interevent Intervals

√ <u>Σ</u> № ≜c ≚c ₩ <b>€</b>	Copy Pre and Post	
B Inst Freq IEI BE	IEI columns to	
3.7702 N/A Not found N/A		
30.842 7.86164 127.2 N/A 108.27 18.7266 53.4 N/A		
16.223 54.6448 Properties 2.0874 108.696	Results1:1 - Results	
71 279 3 8506 Cut Ctrl+X	<u> </u> ∠  ∠ X Y XY Ex 1 A 2 4 w > Σ M Ac ≚c s	
Copy         Ctrl+C           38.674         3.81243         Paste         Ctrl+V           75.404         17.8253         Ctrl+V         Ctrl+V		un b_1
07.662 10.917 Clear Del	Pre Peak     Post Peak     Pre IEI     Post IEI     E     Y     Image: A transmission of the constraint	
03.279 11.0619 Delete Ctrl+D 03.431 9.03343 110.7 N/A	2         38.5418         86.9652         127.2         71.7         st Peak         Pre IEI         Post IEI         E         F           3         38.226         37.8616         44.6134         Not found	G
68 838 12 9534 77 2001 N/A	4 24.0615 18.3423 Column Rename 86.9652 127.2 71.7	
	5         28.3004         14.4603         37.8616         53.4         63.1           6         59.7718         107.519         Mame of column D:         18.3423         000         000	
	7 83.5718 47.0145 Name of column D: Previo 14.4603 Kolmogorov-Smirnov Test	×
	8 69.6571 51.0085 Next 107.519 Preprocessing	
	10 23.0717 37.329 51.0085 6 Ascending sort	
	10         22.011         Gr.325           11         99.0512         21.1148         OK         Cancel         Help         114.19           12         95.6512         13.25         OK         Cancel         Help         114.19           13         95.6512         13.25         OK         Cancel         Help         114.19	
	12 63.3042 13.61 21.1148 © Fixed width: 0.05	
	13.6	
	26.3616	
	100.654 33.8393 C Fixed From: 1 To: 5	
	36.1158	
	37.1242 Select Columns Sheet 18 64.0997	
	44.457 Row Range	
	59.7216         Image: Full column           21.8769         Image: Row         From: Image: Full column	051
	25.6714	851
	Be sure to Append, 67.9402 • Append results to sheet and graph	
	not Replace the 19.5468 C Replace results in sheet and graph	
		Help



#### Together through life sciences.

### **Kolmogorov-Smirnov Test Results**



#### Together through life sciences.



### **Retain the Kolmogorov-Smirnov Results**

🖹 Re	Results1:1 - Results									
$[\underline{\mathcal{M}}]$										
	Source She	Source Colu	n1	n2	K-S Statistic	Probability	Num Bins 1	Num Bins 2		
1	Sheet 18	Pre Peak vs.	334	851	0.07395	0.14510		N/A		
2	Sheet 18	Pre IEI vs. Po	334	851	0.43035	0.00	n Properties			
3							Propercies			
4							Cut (	Etrl+X		
5							Copy (	Etrl+C		
6							Paster (	trl+V		
1 1	<sup>7</sup> Ⅰ ▲ ▶ ਸ਼ & Statistics A Basic Stats & Fit Params & Correlation & Fluctuation & Histog							Del ce <b>,{</b> ∨-		
		 اینا اصحاب					Delete (	Itrl+D		

Copy the results of the two K-S tests to Sheet 17

Results1:1 - Results										
$\boxed{ \mathcal{M} } \boxed{X} \boxed{Y_{+}} \boxed{XY} \boxed{\mathbb{E}_{x}} \boxed{\mathbb{E}_{z}} \boxed{\mathbb{E}_{z}$										
	A B C D E F									
1	Sheet 18 🦵		<u>~</u> 84	851	0.073946	0.145096				
2	Sheet 18	Properties	34	851	0.430346	0				
3		Cut	Ctrl+X							
4		Сору	Ctrl+C							
5		Paste	Ctrl+V							
6		Clear 🕅	Del							
7		Delete	Ctrl+D -							
8	Ļ									
	🕨 🕅 🖌 Shei	et 15 🖌 Sheet 1	16 <mark>)</mark> Sheet 1i	7 🖌 Sheet 18	K Pre K Pos	t/ •				



#### Together through life sciences.

### **Basic Statistics for Peaks and Inter**event Intervals

🖉 Clampfit - Results1												
File Edit View Analyze Format Event Detection Configure Too			Too									
	Create Graph Assign Plots Append to Graph											
Results1:1	Fit Fitting Results		asic Statistic	:s			3					
Image: Market with the second secon	Extract Data Subset Histogram Sort Transpose Average Columns	t IEI ind 71.7	– Measurements Number pe Minimum Maximum Variance Mean		Standard deviation Standard error Coefficient of variation Sum Sum of squares							
3 38 4 24. 5 28.	Column Arithmetic Create Data	63.1 28.7 85.4	Median		Jun of squares							
6 59. 7 83.	Analysis of Variance Basic Statistics	13.7 20	- Destination Op	results to sheet	C Replace results in	sheet						
8         69.           9         92.           10         23.           11         99.           12         85.           13         106           14         32.           15         48.           16         21.           17         36.7           18         422.8	Chi-Square and Mann-Wr&ney F-Test and Student's t-Test Kolmogorov-Smirnov Test Rank Correlations Autocorrelation Cross-correlation Event Analysis 403 36.1158 62.9 692 37.1242 28.2	231 24.9 25.2 62.5 62.9 20.3 60.1 47.9 4.2001 15.2 28.5999	Row Range Full colu Row C Perform Category	elect Columns Sheet 18 All selected columns w Range Full column Row From: 1 To: 1 To: Perform Breakdown Analysis Category column: A: Pre Peak Specify Bins								
	neet 15 <mark>(</mark> Sheet 16 ( Sheet 1 sults1:1 - Results	/ ASheet 18			1	1						
	14 X Y+ XY		A2	»,√  Σ	M <sub>4</sub> ≜c ≚c	1 H						
	Source	#/Cat	Min	Max	Mean	Median	Std. Dev.	Std Err.	CoVar.	Sum	Sum Sq.	
1	Pre Peak:Sheet 18	334	9.61644	247.951	52.2991	39.440	37.6507	2.06015	0.719911		1.38561e+00	
2	Post Peak:Sheet 18	851	8.74508	197.461	48.8093	42.044	31.1591	1.06812	0.638385		2.85264e+00	
3	Pre IEI:Sheet 18	334	0	476.301	89.603	64.300	77.6239	4.24739	0.866309		4.68807e+00	Ļ
4	Post IEI:Sheet 18	851	0	133.399	35.2031	29.200	21.3367	0.731413	0.606104	29957.8	1.44157e+00	╞
5												┝
		1							1			

#### Together through life sciences.

©2012 For research use only. Not for use in diagnostic procedures. Trademarks mentioned herein are property of Molecular Devices, LLC or their respective owners.



Var

1417.57


970.892

6025.47


455.255

### Summary

- Are there differences in the amplitude distribution?
- No.
  - Post-synaptic mechanisms are unaffected.



- Are there differences in the frequency?
- Yes.
  - Pre-synaptic mechanisms are upregulated.





#### Together through life sciences.

## **Clampfit Features Discussed**

- Event Detection > Create Template
- Event Detection > Template Search
- Move Cursors efficiently
- View > Zoom > Between Cursors
- Analyze > Arithmetic
- Format > Column > Rename
- Format > Rename Sheet
- Analyze > Histogram
- Efficiently creating a graph using X-Y pairs
- Editing a graph
- Analyze > Kolmogorov-Smirnov Test
- Analyze > Basic Statistics



#### Together through life sciences.

# **End** Thank You







# Together through life sciences.

www.moleculardevices.com